









COVID-19 Supplemental Oxygen Instructions

Oxygen Administration Device	Image	Protection	
		Type of Mask	Type of Room
Simple face mask 5 LPM		Surgical	Private
Nasal cannula 0.5–6 LPM		Surgical	Private
Oxymask 1–15 LPM		Surgical	Private
Ventimask 24–50%		Surgical	Private
Non-rebreathing mask (filtered)		Surgical	Private
Non-rebreathing mask		Surgical	Private
Airvo / Optiflow	 	N95	Negative pressure